



# GILSON TIPS SELECTION GUIDE

We've created this interactive matrix to help you in selecting the right tip for your instrument and application.

Click on the symbols ● ◆ ■ □ ○ to go directly to the product pages.

FOR PRODUCT DETAILS  
CLICK ON THE SYMBOLS

Purity Grade	Pipette Compatibility	Filter	Packaging Type (Code Description)	Volume Range																			
				0.1-10 µL	0.5-10 µL	0.1-20 µL	1-20 µL	2-30 µL	1-100 µL	10-100 µL	1-200 µL	2-200 µL	10-200µL	20-200 µL	1-300 µL	10-300 µL	100-1000 µL	50-1200 µL	50-1250 µL	0.5-5 mL	1-10 mL		
<b>PIPETMAN DIAMOND TIPS</b> Sterile SAL 10 <sup>-6</sup> * • Free of detectable DNA, RNA, DNase, RNase, ATP, endotoxins (pyrogens) • Trace metal release tested	PIPETMAN	With filter	TIPACK (DFxxxST)	●				●■□○		●■□○				●■□○		●■	●						
			BLISTER REFILL (DFxxxST)	●◆				●◆□○		●◆□○				●◆□○		●◆	●◆						
			BLISTER REFILL Certified 384-well (DSFxxxST)	◆				◆□○						◆□○		◆							
		Without filter**	TIPACK (DxxxST)	●		●							●■□○		●■	●	●						
			TOWERPACK (DxxxST)	●		●										●	●						
			STERILPACK (DxxxST)														●						
AmpliPur Expert Tips Sterile SAL 10 <sup>-6</sup> * Free of Human DNA, RNA, DNase, RNase	Most leading brands	With filter	TIPACK		●		●									●	●						
			PIPETMAN EXPERT TIPS Sterilized SAL 10 <sup>-3</sup> * • Free of human and bacterial DNA, DNase, RNase, non pyrogenic, no inhibition of real time PCR Available in selected countries only	Most leading brands	Without filter**	TIPACK (EFxxxST)	●		●			●					●						●
TIPACK (ExxxST)	●		●							●					●						●		
TIPACK - Gel loading flat or round tips (ExxxST)**											●												
PIPETMAN DIAMOND TIPS Autoclavable** • Manufactured in clean room environment • Trace metal release tested	PIPETMAN	Without filter	TOWERPACK (Dxxx)	●		●									●	●							
			TIPACK (Dxxx)	●		●							●■□○		●■□	●	●			●	●		
			BLISTER REFILL (Dxxx)			●◆								●◆□○		●◆□	●◆	●◆					
			BLISTER REFILL Certified 384-well (DSxxx)			◆								◆□○			◆						
PIPETMAN DIAMOND TIPS Autoclavable** • Trace metal release tested	PIPETMAN	Without filter	EASYPACK (Dxxx)	●		●								●■□	●	●			●	●			
			RELOAD PACK (Dxxx)			●							●■□○			●							
			ECOPACK (Dxxx)	●		●							●■□			●■□	●	●					
PIPETMAN EXPERT TIPS Autoclavable** • Standard Available in selected countries only	Most leading brands	Without filter	TOWERPACK (Exxx)**			●									●					●			
			TIPACK (Exxx)	●		●									●						●		
			TIPACK - Gel loading flat or round tips (Exxx)**	●													●						
			EASYPACK (Exxx)			●																	

## Tips Legend

- Compatible with PIPETMAN\*
- ◆ Compatible with PIPETMAX\* 268
- Compatible with PLATEMASTER\* P20
- Compatible with PLATEMASTER\* P220
- Compatible with PIPETMAN\* M96 P20
- Compatible with PIPETMAN\* M96 P200



PIPETMAN® DIAMOND Tips are designed to fit PIPETMAN® perfectly and guarantee maximum accuracy and precision in conjunction with your Gilson pipette.

PIPETMAN® EXPERT Tips can be used with most leading pipette brands.

\*Sterility Assurance Level SAL (ISO 11137): quantitative value to assure sterility. SAL 10<sup>-6</sup>: probability of 1 in 1.000.000 of finding a non-sterile unit.

\*\*Autoclavable: tips are fully autoclavable except gel loading tips. Autoclaving does not destroy traces of DNA and RNases. Sterile not filtered PIPETMAN DIAMOND Tips and Sterilized not filtered PIPETMAN EXPERT tips are autoclavable.

TOWERPACK Tips are autoclavable, but the plastic wrap and the tray of the TOWERPACK are not. We recommend to place the racks in the tip boxes before autoclaving.

Empty reload boxes and racks are available for RELOAD, BLISTER REFILL, TOWERPACK, ECOPACK, EASYPACK packaging options.

## WIDE RANGE OF TIP PACKAGING

### PIPETMAN® DIAMOND TIPS

#### TIPACK™

The design is more compact, convenient, ergonomic and eco-friendly

Patented ROCKY RACK system and hinged 180° angled lid makes it extremely easy to fit tips on multichannel pipettes

Fully autoclavable

Available as sterilized and filter tips

[VIEW ALL](#)



#### TOWERPACK

A quick, easy, and affordable rack refill system

Uses 60% less space on your bench and 70% less plastic than regular rack boxes

The whole range of TOWERPACK fits the amazingly easy-to-use universal reload box

One-hand system operation

Fully autoclavable

Available sterilized

[VIEW ALL](#)



#### BLISTER REFILL

Individually wrapped racks of tips that are convenient, clean, and eco-friendly

Ready to use, just peel the lid and insert into the empty box

Uses 45% less plastic than regular rack box

Fully autoclavable

Available as sterilized and filter tips

[VIEW ALL](#)



#### EASYPACK AND ECOPACK

Purchase in bulk and save money by simply reloading these tips into empty racks

Bags are sealed with rapid-sealing strip for easy access and to prevent contamination

Minimal waste

Fully autoclavable

EASYPACK: 5 bags x 200 tips

ECOPACK: 10 bags x 1000 tips

[VIEW ALL](#)



#### STERILPACK

Individually wrapped sterilized tips

Complete protection from contaminants up to the moment of use

Wrapped in medical-grade quality plastic

Sterilization guaranteed for 5 years

400 individually wrapped sterilized tips

[VIEW ALL](#)



#### RELOAD PACK

Quickest way to reload empty box

Uses 70% less plastic, recyclable packaging is eco-friendly

Easy to store, takes up less bench space

Affordable and convenient

Fully autoclavable

[VIEW ALL](#)



### PIPETMAN® EXPERT TIPS

#### TIPACK™

Filter, autoclavable, sterilized

Uniquely designed TIPACK features a dual-purpose lid

Hinged to easily open and close while using tips

Pull it out of the hinges to set on top or discard for high-throughput tip usage

[VIEW ALL](#)



#### MEGAPACK

Autoclavable, sterilized

Uniquely designed TIPACK features a dual-purpose lid

Sterilized MEGAPACKS give you the choice to expose 96 tips at a time or 192

Ideal for high throughput pipetting

[VIEW ALL](#)



#### EASYPACK

Autoclavable

Purchase in bulk and save money by simply reloading these tips into empty racks

Resealable bags for easy access

EASYPACK: 5 bags x 200 tips

[VIEW ALL](#)



#### TOWERPACK

A quick, easy and affordable rack refill system

Uses 60% less space on your bench and 70% less plastic than regular rack box

One-hand system operation

[VIEW ALL](#)



### AMPLIPUR® EXPERT TIPS

#### AMPLIPUR® EXPERT TIPS

Filter tips, DNA, RNA, DNase and RNAase free

Best value for money

Ready-to-use and optimized for use with multichannel

Pipettes with the rack lid opening position at 180°

[VIEW ALL](#)

