



VARIABLE VOLUME MECHANICAL PIPETTES



I Ultimate Performance

Built for traceability and comfort. PIPETMAN L is designed to increase daily productivity, giving you the ability to pipette longer without fatigue.

I Secure Pipetting

Reduce pipetting errors with volume-locking. Avoid accidental volume changes during pipetting cycles, ensuring your accuracy and precision.

I Optimized Comfort

Lightweight and light pipetting forces. Right-handed or left-handed, no problem. Adjust the location of the tip ejector to the most comfortable position for you.





Classic metal or plastic tip ejector

Choose between a stainless steel tip ejector for more robustness or a plastic tip ejector to avoid corrosion.



Choose your model

Choose your channel number and volume range. Each ensures consistent volume levels, allowing you to aspirate and dispense your samples with confidence.

PIPETMAN L Single Channel — Variable Volume Models				
Model	PIPETMAN DIAMOND Tips		Part Number	Volume (µL)
P2L	D10 DL10	DF10 DFL10	FA10001P FA10001M	0.2 - 2
P10L*	D10 DL10	DF10 DFL10	FA10002P FA10002M	1 - 10
P20L	D200	DF30	FA10003P FA10003M	2 - 20
P100L	D200	DF100	FA10004P FA10004M	10 - 100
P200L	D200 D300	DF200 DF300	FA10005P FA10005M	20 - 200
P1000L	D1000 D1200	DF1000 DF1200	FA10006P FA10006M	100 - 1000
P5000L	D5000		FA10007**	500 - 5000
P10mLL	D10mL		FA10008**	1000 - 10000
PIPETMAN L Multichannel				
P8x10L	D10 DL10 DF10 DFL10	FA10013		0.5 - 10
P12x10L		FA10014		
P8x20L	DL10 DFL10 DF30 D200	FA10009		2 - 20
P12x20L		FA10010		
P8x200L	D200 D300 DF200 DF300	FA10011		20 - 200
P12x200L		FA10012		
P8x300L	D200 D300 DF200 DF300	FA10015		20 - 300
P12x300L		FA10016		

For a pipette with a plastic tip ejector, the part number ends with the letter P, For a pipette with a stainless steel tip ejector, the part number ends with the letter M.

PIPETMAN L specifications are downloadable from www.gilson.com/Specifications/PipetmanL

PIPETMAN L single channel models are available in kits.

For even more reproducibility, PIPETMAN L is also available in a fixed version.

www.gilson.com/ContactUs



LT800565 © 09/2017 Gilson S.A.S.
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

