


# pipetman<sup>®</sup>

L



## INCREASED COMFORT AND PRODUCTIVITY

Air-Displacement Pipettes



**Designed for performance and comfort, PIPETMAN<sup>®</sup> L offers a comprehensive range of models, adapted to all your needs.**

### Convenient Use

Lightweight and balanced, PIPETMAN L combines an ergonomic design with minimized pipetting forces, giving you the ability to pipette longer without fatigue.

Right-handed or left-handed? No problem. Adjust the location of the tip ejector to the most comfortable position for you.

### Secure Pipetting

Reduce pipetting errors and avoid accidental volume changes during pipetting cycles with a true volume-locking mechanism for improved accuracy and precision.

PIPETMAN L variable volume models are all fully autoclavable without disassembly for maximum convenience and protection from cross-contamination.

### Extended Range

Available as a fixed or variable volume, the PIPETMAN L line covers a full volume range from 0.2  $\mu\text{L}$  to 10 mL for single-channel models, and from 0.5  $\mu\text{L}$  to 1200  $\mu\text{L}$  for multichannel models. Precise, accurate, and reliable, they can be intensively used with confidence in many applications.

Four PIPETMAN L multichannel models (product description ending with VR), covering a volume range from 20  $\mu\text{L}$  to 300  $\mu\text{L}$ , are equipped with V-rings at the bottom of the tip holders for a proper fit with most tip brands.



Adjustable Tip Ejector



## Choose your Tip Ejector

For single channel models, choose between a model with a stainless-steel tip ejector or a plastic ejector depending on your needs and applications.\*

Large volume models are equipped with a plastic tip ejector to avoid any direct contact with the tip.



## Choose your Model

Choose your channel number and volume range, up to 1200 µL for various protocols.

Four PIPETMAN L multichannel models equipped with V-rings (VR) offer an improved compatibility with tips from other brands.



## Choose your Pipette Stand

Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.



PIPETMAN L single channel models are available in kits, offering the essential pipetting tools to start an experiment right away.



For routine testing applications, PIPETMAN L is also available as a fixed volume version, covering a volume range from 1 µL to 5 mL.



PIPETMAN DIAMOND® Tips are designed to fit your PIPETMAN perfectly and guarantee maximum accuracy and precision.

For ordering information, visit [www.gilson.com](http://www.gilson.com)

### PIPETMAN L Single Channel – Variable Volume Models

| Model  | PIPETMAN® DIAMOND Tips |                      | Part Number | Volume (µL)  |
|--------|------------------------|----------------------|-------------|--------------|
| P2L    | D10<br>DL10            | DF10ST<br>DFL10ST    | FA10001P    | 0.2 - 2      |
|        |                        |                      | FA10001M    |              |
| P10L   | D10<br>DL10            | DF10ST<br>DFL10T     | FA10002P    | 0.5 - 10     |
|        |                        |                      | FA10002M    |              |
| P20L   | D200                   | DF30ST               | FA10003P    | 2 - 20       |
|        |                        |                      | FA10003M    |              |
| P100L  | D200                   | DF100ST              | FA10004P    | 10 - 100     |
|        |                        |                      | FA10004M    |              |
| P200L  | D200<br>D300           | DF200ST<br>DF300ST   | FA10005P    | 20 - 200     |
|        |                        |                      | FA10005M    |              |
| P1000L | D1000<br>D1200         | DF1000ST<br>DF1200ST | FA10006P    | 100 - 1000   |
|        |                        |                      | FA10006M    |              |
| P5000L | D5000                  |                      | FA10007     | 500 - 5000   |
| P10mLL | D10mL                  |                      | FA10008     | 1000 - 10000 |

\*For a pipette with a plastic tip ejector, the part number ends with the letter P.

For a pipette with a stainless steel tip ejector, the part number ends with the letter M, from P2L to P1000L.

### PIPETMAN L Multichannel Models

| Model      | PIPETMAN® DIAMOND Tips |                    | Part Number | Volume (µL) |
|------------|------------------------|--------------------|-------------|-------------|
| P8x10L     | D10<br>DL10            | DF10ST<br>DFL10ST  | FA10013     | 0.5 - 10    |
| P12x10L    |                        |                    | FA10014     |             |
| P8x20L     | DL10<br>D200           | DFL10ST<br>DF30ST  | FA10009     | 2 - 20      |
| P12x20L    |                        |                    | FA10010     |             |
| P8x200L    | D200<br>D300           | DF200ST<br>DF300ST | FA10011     | 20 - 200    |
| P12x200L   |                        |                    | FA10012     |             |
| P8x200LVR  | D200<br>D300           | DF200ST<br>DF300ST | FA10035     | 20 - 200    |
| P12x200LVR |                        |                    | FA10036     |             |
| P8x300L    | D200<br>D300           | DF200ST<br>DF300ST | FA10015     | 20 - 300    |
| P12x300L   |                        |                    | FA10016     |             |
| P8x300LVR  | D200<br>D300           | DF200ST<br>DF300ST | FA10037     | 20 - 300    |
| P12x300LVR |                        |                    | FA10038     |             |
| P8x1200L   | D1200                  | DF1200ST           | FA10039     | 100 - 1200  |
| P12x1200L  |                        |                    | FA10040     |             |

### Pipette Stands

| Description                          | Part Number |
|--------------------------------------|-------------|
| CARROUSEL Pipette Stand (7 pipettes) | F161401     |
| TRIO Pipette Stand (3 pipettes)      | F161405     |
| Universal Multichannel Stand         | F161417     |
| SINGLE Pipette Holder                | F161406     |

[gilson.com/contactus](http://gilson.com/contactus)

**ISO 9001**  
CERTIFIED

LT800606/A | ©2021 Gilson S.A.S.  
All rights reserved. All trademarks  
are the property of Gilson, Inc.  
and its subsidiaries.

**GILSON®**